Advanced Bionics AG
Laubisrütistrasse 28
8712 Stäfa, Switzerland
+41 58 928 78 00

Manufactured by:
Advanced Bionics LLC
California, U.S.A.
+1 661 362 1400

AdvancedBionics.com

029-M486-85 Rev A

© 2015 Advanced Bionics AG. All Rights Reserved.

MRI Antenna Coil Cover

CI-7521

instructions for use

English, Français, Español, Deutsch, Nederlands, Português

Labeling symbols and their meanings:

 REF	Model number		Date of manufacture
 LOT	Lot Number		Manufacturer
	Do Not Reuse		See Instructions for Use
	MR Conditional: Specific circumstances are required for the patient to undergo an MRI procedure with this device in place. See complete Instructions for Use for details		

MRI Warnings

MRI is contraindicated except under the circumstances described below. Do not allow patients with a HiRes 90K or HiRes 90K Advantage cochlear implant to be in the area of an MRI scanner unless the following conditions have been met:

- The bandaging protocol recommended by Advanced Bionics is followed when the patient undergoes an MRI procedure with the magnet left in place or,*
- the internal magnet is surgically removed and possibly replaced with the Magnet Insert Dummy before the patient undergoes an MRI procedure.*

- *The external sound processor and headpiece are removed before entering a room where an MRI scanner is located.*
- *MRI preparations must be conducted outside of the MRI room.*
- *Verify that the implant, or both implants if bilaterally implanted, are compatible for conducting an MRI before proceeding. Failure to do so can lead to device movement, device damage, magnet movement, patient discomfort, or trauma and pain to the patient.*
- *MRI procedures are contraindicated for CLARION (C1 and CII) Cochlear Implant recipients.*

Intended Use

The MRI Antenna Coil Cover is intended to be used along with a bandaging protocol to allow the patient to undergo an MRI procedure with the magnet left in place. Instructions for use must be followed, including the recommended bandaging protocol.

MRI Antenna Coil Cover

Compression Dressing prior to MRI Procedure with HiRes 90K Users

This document will provide general guidance for applying appropriate compression dressing over the site of the implanted device (and magnet) in order to perform an MRI procedure.

Required Materials


- Hurt-Free Wrap – Hospital Grade
Latex Free 5 cm width
Johnson and Johnson (or equivalent)
- Flexible tape
Latex Free
2 – 2.5 cm width
3M – Nexcare (or equivalent)
- Patient Headpiece
- MRI Antenna Coil Cover
CI-7521
Advanced Bionics
- Marking Pen

Procedure


1. Place patient in sitting position to allow access to the implant site.

NOTE: *If the patient is bilaterally implanted, prepare to locate both implant devices using the following steps. Both implant locations must have compression dressing placed over their magnets.*


2. Obtain the patient headpiece (with cable removed) and place over the implant site. Allow magnetic attraction to occur.

3. Once in place, hold the headpiece while drawing a line around it.

4. Obtain the MRI Antenna Coil Cover and place it over the drawn location. Once centered, place flexible tape over the site to maintain the location.

5. Wrap the head using Hurt-Free Wrap sufficiently tight to ensure stabilization of the MRI Antenna Coil Cover.

After the MRI

Once the MRI procedure has been completed, gently remove the Hurt- Free Wrap, tape, and the MRI Antenna Coil Cover.

Discard all components including the MRI Antenna Coil Cover.

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