
Labeling symbols and their meanings:

 0123	European Community Mark of Conformity. Authorized to affix CE Marking in 2013	
	Model number	 Date of manufacture
	Serial Number	 Manufacturer
	Type of Protection: B	 Store at temperatures between -20°C (-4°F) and $+55^{\circ}\text{C}$ (131°F)
	Fragile	 Do not get wet
	Lot Number	 See Instructions for Use
	Suitable for atmospheric range between 70kPa and 106kPa, which is equivalent to 3000m above sea level to 380m below sea level.	
	Suitable for exposure to relative humidity between 0-95%	
	Dispose of in accordance with applicable national and local regulations	


Warnings

- **CHOKING HAZARD** –The Naída CI contains small parts that pose a hazard of inhalation or choking
- It is recommended to use the T-Mic 2 only when the T-Mic 2 Cover is properly seated.

Naída CI T-Mic 2 Microphone and Earhook

The Naída CI T-Mic 2 and earhook are intended for use with the Naída CI sound processor. The T-Mic 2 is a custom in-the-ear microphone for hearing in noisy environments and gaining effortless access to cell phones, consumer audio electronics, and assistive listening devices while the earhook serves to hold the processor on the ear when the built-in microphone is utilized.

Connecting the T-Mic 2 or Earhook

Beginning with the T-Mic 2 or earhook disconnected, align the T-Mic 2 or earhook flush with the processor. Hold the T-Mic 2 or earhook firmly at the top and push towards the processor until it meets flush with the processor.

Using the provided T-Mic 2 or earhook tool, place one of the pins in the small hole on the side of the tool. This will help to steady the pin as you line it up with the pin hole on the T-Mic 2 or earhook.

Line up the pin with the hole in the T-Mic 2 or earhook and gently push it through the hole. The pin will move from one side of the processor, through the T-Mic 2 or earhook to the other side of the processor, holding the T-Mic 2 or earhook in place. The pin should not come out the other side.


If a small portion of the pin extends beyond the processor, use the side of the tool to push the pin inwards so it lies flush with the processor.


Pull the tool away from the processor, leaving the pin inside the Naída CI.


Removing the T-Mic 2 or Earhook

Remove the power source. Align the pointed end of the 2 tool with the hole on the processor that contains the pin holding the T-Mic 2 or earhook on the Naída CI.

Gently press the tool in the hole, pushing the pin out the other side.

The pin does not need to be removed from the processor. As long as the tool was pushed in flush on one side of the processor, the pin should be displaced enough to adequately remove the T-Mic 2 or earhook.


Remove the tool from the Naída CI. Gently grasp the T-Mic 2 or earhook at the top near where it connects with the processor and pull forward until it is separated from the processor.

Attaching the T-Mic 2 Cover*

Hold the T-Mic 2 Cover. In the other hand, hold the T-Mic 2 at the end that houses the microphone. Align the bottom of the T-Mic 2 over the T-Mic 2 Cover.

Gently turn the T-Mic 2 into the T-Mic 2 Cover. Screw until they meet together flush.


Once you feel resistance, stop turning.

Removing the T-Mic 2 Cover


Grasp the end of the T-Mic 2 where the T-Mic 2 Cover is located, in one hand. Firmly hold the T-Mic 2 microphone portion in the other hand.

Note: The T-Mic 2 Cover should be replaced at regular intervals to maintain sound quality.

Turn the hand holding the T-Mic 2 Cover to slowly unscrew the T-Mic 2 Cover from the T-Mic 2.

Continue turning until the T-Mic 2 Cover separates from the T-Mic 2.

NOTE: Do not attempt to service or modify the Naida CI or its accessories. Doing so may compromise system performance and will void the manufacturer's warranty. Products should be serviced only at Advanced Bionics.

Recommended Operating and Storage Temperature Ranges

Condition	Minimum	Maximum
Operating Temperature	0°C (32°F)	45°C (115°F)
Storage Temperature	-20°C (-4°F)	55°C (131°F)

Naída CI External Equipment and Accessories

The products listed below are compatible

Description	Model Number
Naída CI Q90 Sound Processor* [†]	CI-5280-120, CI-5280-130, CI-5280-140, CI-5280-150, CI-5280-160, CI-5280-170, CI-5280-190, CI-5280-110, CI-5280-220, CI-5280-230, CI-5280-240
Naída CI Q70 Sound Processor* [†]	CI-5245-110, CI-5245-120, CI-5245-130, CI-5245-150, CI-5245-160, CI-5245-170, CI-5245-140, CI-5245-190, CI-5245-220
Naída CI Q30 Sound Processor* [†]	CI-5260-120, CI-5260-140, CI-5260-150
Universal Headpiece (UHP)*	CI-5305
Concave UHP*	CI-5307
T-Mic 2 Microphone*	CI-5835-100, CI-5835-200, CI-5835-300
Earhook*	CI-5710-100, CI-5710-200


*Applied part per IEC 60601-1

[†]Only available in markets where regulatory approval has been received and the product has been made

commercially available. Contact Advanced Bionics for more information.

NOTE: model number extensions denote variants related to product color, size, and/or length.

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