



How the Ear Works

For sound to be heard, the ear needs to convey the message to the brain. The ear consists of four main parts: outer ear, middle ear, inner ear, and the auditory nerve. Each of these parts plays a key role in sound transmission to the brain.

The Outer Ear

The outer ear (or pinna) captures sound and then directs the sound down the ear canal to the eardrum.

The Middle Ear


Sound vibrates the eardrum and is transmitted to the middle ear. The middle ear consists of three bones, the malleus (hammer), incus (anvil), and stapes (stirrup) all of which help convey sound to the inner ear.

The Inner Ear

The inner ear (or cochlea) consists of sensory receptor cells. These sensory cells have distinct functions in the transmission of sound. It is the unique properties of the inner ear that are responsible for transcribing pitch and loudness information into an electrical code that the nerve can then transmit to the brain.

The Auditory Nerve

The auditory nerve (or hearing nerve) is responsible for transmitting the sound information from the inner ear to the higher processing centers of the brain. The brain is responsible for interpreting sound.



What is “Normal” Hearing?

When the ear mechanisms work properly, the ear can detect a range of pitches or frequencies over a wide range of loudness or intensity. Frequency is measured in Hertz (Hz) and intensity in decibels (dB). Normal-hearing persons can detect very soft sounds of 20 dB to very loud sounds of 120 dB, over a wide range of frequencies.

An audiogram plots a person’s threshold—the softest level the person can hear as a function of frequency and intensity. People with thresholds of 20 dB or below are considered to have normal hearing.



What is Hearing Loss?

Hearing loss refers to a decrease in a person's sensitivity to sound. Alternatively, hearing loss is an increase in a person's thresholds to sound. Typically, hearing loss involves damage to the sensory cells of the inner ear, referred to as a sensorineural hearing loss. Hearing losses will range in degree from mild to profound, depending on the extent of sensory cell loss or damage.

Mild Hearing Loss is designated by hearing sensitivity in the 25–40 dB range. Moderate Hearing Loss is designated by hearing sensitivity in the 40–70 dB range. Severe Hearing Loss is designated by hearing sensitivity in the 70–90 dB range. Profound Hearing Loss is designated by the 90 dB and greater range. Often, persons with profound hearing loss are referred to as deaf.

Damage to the sensory or neural structures result in a permanent hearing loss. There is no medical cure for sensorineural hearing loss.


What is a Hearing Aid?

For persons with hearing losses, the most commonly recommended treatment to improve hearing sensitivity is a hearing aid. A hearing aid amplifies sound, or more simply, makes sound louder. By amplifying sound, any remaining, healthy sensory cells in the ear are stimulated to transmit sound information to the auditory nerve and brain.

What is a Cochlear Implant?

For patients with severe-to-profound sensorineural hearing losses, hearing aids are often not able to improve speech understanding. While very powerful hearing aids exist that increase intensity to very high levels, they may not improve hearing if there are too few or no remaining sensory cells to stimulate. In these cases, the sound information cannot be properly delivered to the brain. For these patients, a cochlear implant may be recommended. Cochlear implants attempt to restore hearing in people with severe-to-profound hearing loss by bypassing the damaged structures of the inner ear and delivering electrical stimulation directly to the auditory nerve. Cochlear implants are prosthetic devices with internal components that are surgically placed and external components that require fitting and programming. Cochlear implants are considered the only medical treatment for severe-to-profound hearing loss.

Cochlear implants attempt to restore hearing to people with severe-to-profound hearing loss by bypassing the damaged structures of the inner ear and delivering electrical stimulation directly to the auditory nerve.





How is a Cochlear Implant Different from a Hearing Aid?

Hearing aids acoustically amplify sound and rely on the responsiveness of healthy inner ear sensory cells to receive that sound and send the message to the brain. In patients with severe-to-profound sensorineural hearing loss, the inner ear sensory cells are lost or damaged to the degree that a hearing aid is no longer beneficial for understanding speech (i.e., there are not enough healthy sensory cells to relay the message to the brain). A cochlear implant, however, bypasses the absent or damaged sensory cells and stimulates the hearing nerve directly by converting acoustic sound input into an electrical pattern that is recognizable to the hearing nerve for the transmission of sound information to the brain.

Who is a Candidate for a Cochlear Implant?

Adults and children 12 months of age and older with severe-to-profound sensorineural hearing loss may be candidates for a cochlear implant. Typically, candidates will have been fitted with hearing aids but received minimal benefit. Benefit with hearing aids is determined by standard diagnostic tests administered by audiologists, as well as parent and teacher/therapist questionnaires, which measure a child's response to sound and development of speech with hearing aids.

A team approach is most often taken to determine if a child is a candidate. A surgeon rules out medical contraindication with a physical examination, an MRI or a CT scan, and medical review of records. An audiologist evaluates the degree of hearing loss, functional benefit with hearing aids, and the potential benefit to be gained from a cochlear implant. A speech-language pathologist evaluates communication development, including speech and language skills. The psychologist assesses cognitive, sensory, and motor development as well as parental expectations and motivations. Educators provide information regarding the child's academic progress. The key member of the team is the parent(s)/caregiver(s) whose observations are essential in evaluating the child's overall development and communicative function in the home and outside the educational setting.

The evaluation allows the team to share experiences and knowledge with the family so that appropriate goals and expectations regarding outcome and rehabilitation can be made. If the team determines the child is a candidate, each member of this team will also play a crucial role in the child's success with the device.

What is the Procedure for Getting a Cochlear Implant?

Preparing the Child for Surgery

Once a child is determined to be a candidate for a cochlear implant, it is helpful to prepare the child for the surgical placement of the internal portion of the device. It is recommended that the child be familiarized with the surgical setup—clothing (caps, gowns, and masks) and the anesthesia mask. Many hospitals or implant centers have materials and programs to assist families and educators in familiarizing the child with this process. Advanced Bionics sends either the Bionic Buddy or Melody, our monkey mascots, with each child's implant system. Bionic Buddy and Melody wear a sound processor(s) and may be useful in explaining the device to the child.

Surgical Placement of Internal Components

The surgical placement of the internal portion of the device is performed under general anesthesia and takes approximately two hours. Typically, the surgery is done on an outpatient basis in which the child arrives in the morning and goes home that same day. Occasionally, an overnight stay may be necessary.



Fitting the External Components

Following a three-to-six-week recovery period, the child will return to the audiologist for fitting of the external components. This first fitting is often referred to as the initial stimulation or hook-up. An audiologist programs the device during the fitting process. Programming involves setting specific parameters of stimulation for the recipient, particularly the levels perceived as soft and comfortable. Programming children is different than programming adults. Children typically do not have the language or experience with sound to report when stimulation is perceived as soft or loud. In addition, children do not have the attention span to sit for an hour of programming. Audiologists who specialize in pediatrics use their experience to fit the device for this more challenging population.

Programming requires frequent adjustments or fine-tuning as the child adapts to the device. While some children adjust very quickly, others require several weeks or months. How much sound a child understands depends on several factors, including age at implantation, length of deafness, previous experience with sound, and access to aural rehabilitation and therapy services.

Is Aural Rehabilitation Required for Cochlear Implant Recipients?

Aural rehabilitation or therapy is a key factor for successful use of a cochlear implant. The implant is not a cure for hearing loss. A new user will need to practice listening with the device to reach maximum performance. For children, rehabilitation is crucial for the development of oral language and speech skills. Parents of children with hearing loss will have several options in education and therapy methodologies. Some of these options use manual forms of communication, such as sign language, while others emphasize oral communication and/or lipreading. Regardless of which educational approach, some form of aural rehabilitation or listening practice is needed to maximize performance. For further information on how to maximize a child's success with a cochlear implant, visit AdvancedBionics.com and click on 'Support Library' or 'Education' to view recorded webinars. Additionally, we encourage you to visit The Listening Room™ at TheListeningRoom.com, an AB online rehabilitation resource for you, your child and his or her educational team.



TheListeningRoom.com • HearingJourney.com

Advanced Bionics
AdvancedBionics.com

For questions or additional information:
Toll Free Phone: 1-877-829-0026
TTY: 1-800-678-3575
CustomerService@AdvancedBionics.com
ToolsForSchoools@AdvancedBionics.com
Monday through Friday, 5am to 5pm PST