

## SoundWave Database Review Tool Download Page

This page includes information on the download, installation, use, and interpretation of the SoundWave Database Review Tool. This tool is intended for use by professionals actively engaged in the ongoing care of Advanced Bionics patients with a HiRes Ultra or HiRes Ultra 3D cochlear implant.

### About the software

This software is compatible with Windows 7 and 10.

The software must be installed on a computer with an active installation of SoundWave 3.2

### Download

The software can be downloaded by clicking this [\[LINK\]](#).

### Installation

The downloaded file is an executable (.exe) file.

The file should be copied to the computer you wish to use for analysis.

The computer must have a Windows operating system (Windows 7 or 10).

The computer must have an active installation of SoundWave 3.2.

### Use

Double-click the icon for the executable file to open the user interface

The software will automatically connect to the database used by SoundWave 3.2

- The connected database is displayed in the interface
- The connected database can be changed by these steps
  - o Close the SoundWave Database Review Tool
  - o Open SoundWave 3.2
  - o Establish a connection to the desired database
  - o Close SoundWave 3.2
  - o Open the SoundWave Database Review Tool
    - The desired database should be connected automatically

Enter the name of your center (required)

Enter the user name (required)

Indicate if patient name should be included in the report (checkbox)

If available, enter a RADIUS Case ID (optional)

Click "Start Analysis" to review files in the connected database

The progress bar will update automatically

Once complete, two reports will be generated automatically

### Interpretation

The detailed report will include a summary of findings and a full Device List

The summary report will only include the Summary Section

The results are based on a custom algorithm used to identify specific impedance patterns associated with a Field Action for HiRes Ultra and HiRes Ultra 3D devices. The report will list devices with high and moderate matches to the search criteria.

#### Summary Section

- Total devices analyzed
- Number of "High Match" devices
- Number of "Moderate Match" devices
- Number of unanalyzed devices due to "not enough measurement data"
- Number of devices where there was "No Problem Detected"

#### Device List

Each device in the database will be included in the appropriate section, based on the analysis.

The section includes patient name (if selected when running the report) and device identifiers. Intellilink ID and Implant Side are extracted from the SoundWave Database. The Serial Number will be available for all V1 devices, but may not be displayed for later versions.

Devices with a “High” or “Moderate” match to known impedance patterns should be reviewed in detail. Specific device concerns should be discussed with your local Advanced Bionics Representative.