

25<sup>th</sup> September, 2017

Dear Hearing Healthcare Professional,

Advanced Bionics is excited to announce the availability of the next version of our fitting software. With a single step, SoundWave™ 3.1 will enable compatibility of Naída CI Q70 and Q90 sound processors with Phonak Naída™ Link CROS, offering your unilateral recipients A Simple Way to Hear More. SoundWave 3.1 also delivers several workflow and efficiency improvements, allowing you to maintain focus where it matters most: on your patients.

What's new in SoundWave 3.1?

- Initialization of the Naída CI Q70 and Q90 enables compatibility with the **Naída Link CROS** in a single step. No fitting of the CROS device is required!
- Automatic management of microphone and auxiliary signals enables the use of **RogerReady** for any Naída program, providing easy access to the powerful technology of the Roger system without the need for program changes.
- Creation of **progressive programs** in your preferred increments is now possible without adding time or complexity to your fitting session.
- New **default preferences** allow the professional to set up default program parameters based on personalized fitting preferences for various age groups.

As part of our commitment to improve the usability and reliability of our CI system, several hardware enhancements for the Naída CI system will also become available throughout the Summer and Fall 2017 seasons. Improved UHP cables will be available in a more streamlined set of cable lengths, making room for the introduction of our new gray cable. A new earhook with grip will also be offered to clinicians who require additional retention when using a standard earmold configuration on the Naída CI. We look forward to letting you know about ongoing system improvements which will be announced once available.

Your regional Advanced Bionics team will provide additional information regarding support and training for SoundWave 3.1, Naída Link CROS, and Naída CI hardware enhancements.

Kind Regards,



Erin Castioni, Au.D.  
Product Manager – Fitting Software