



For more information on hearing your best with AB, visit

AdvancedBionics.com

*Pending regulatory approval in the United States

**Not approved for pediatric use in the United States

028-M064-03

©2011 Advanced Bionics, LLC and its affiliates. All rights reserved.



NEPTUNE™

by Advanced Bionics

The industry's first freestyle™ design

- Hear and communicate in the pool
- Bond during bath time
- Chic, modern freestyle design inspired by consumer electronics
- Compatible with the industry's most advanced hearing technology, including ClearVoice™*, HiRes Fidelity 120™**, and AutoSound™
- Removable controls for waterproof and childproof wearing configuration
- Wide range of fashion-forward colors for mixing and matching
- Two headpiece choices: AquaMic or Universal Headpiece
- Power it with one off-the-shelf AAA battery, disposable or rechargeable
- Convenient wearing accessories: swim cap, armband, headband, clip, pouch, lanyard



UHP AquaMic



Colors	Neptune	UHP	AquaMic
Onyx	●	●	●
Gray		●	●
Dark Sienna		●	
Dark Brown		●	●
Beige Metallic		●	
Beige			●
Sand Castle	●	●	
Ocean Blue	●	●	
Sky Blue	●	●	
Tropical Sea	●	●	
Lime Green	●	●	
Wild Violet	●	●	
Ruby Red	●	●	
Raspberry		●	
Princess Pink	●	●	
Alpine White	●	●	●

colors shown here may vary from final product

NEPTUNE

by Advanced Bionics



wear it your way, every day!

IP 68[†] strength in numbers

AB has the industry's highest Ingress Protection rating!



6

Protection against solid foreign bodies

Rating: 0 indicates no protection and 6 indicates dust tight. AB offers maximum protection against solid particles.



8

Protection against various forms of moisture

Rating: 0 indicates no protection and 8 indicates protection against complete, continuous submersion in water. Only AB offers the world's first processor and microphone intended for repeated submersions.

[†]Complies with IEC 60529 IP68 rating.